



T CELLS – UNDER-APPRECIATED IN COVID-19


RECOGNITION FROM ROUTINE & RESEARCH WITH IN VITRO ASSAYS

Dr. Lester Thoo
ADR-AC GmbH, Bern

12.05.2022 ADR-AC Symposium, Bern


Adaptive immunity in Covid-19

Cellular




Current
diagnostic
tests

Humoral



Cellular immunity test

Cellular



SARS-CoV-2 structure not just Spike!

Large single-stranded
RNA virus (30 kb)

16 non-structural proteins

**3 abundant
structural proteins**

Coronavirus Structure and Protein Visualization






FIGURE 1: Coronaviruses belong in the family *Coronaviridae* and can cause disease in mammals and birds. The coronavirus spike (S) protein mediates membrane fusion by binding to cellular receptors.

Assay Protocol



Flow cytometry
(cell activation markers)

Helper
T cells

Cytotoxic
T cells

-	1	RPMI
-	2	RPMI + DMSO
+	3	Phytohaemagglutinin (PHA)
T	4	Spike 1 protein
T	5	Spike 2 protein
T	6	Nucleocapsid protein
T	7	Membrane peptides mix
T	8	Spike peptides mix
T	9	Nucleocapsid peptides mix

Use in Routine Diagnostics

Established in 1871

Swiss Medical Weekly

Formerly: Schweizerische Medizinische Wochenschrift

An open access, online journal • www.smw.ch

Original article | Published 13 September 2021 | doi:10.4414/SMW.2021.w30005

Cite this as: Swiss Med Wkly. 2021;151:w30005

Highly specific and reliable in vitro diagnostic analysis of memory T and B lymphocytes in a Swiss cohort of COVID-19 patients

Lester Thoo^a, Pierre I. Gumowski^{bc}, Kevin Kammermann^a, Swelia Nussli^a, Benno Grabscheid^a, Oliver Hausmann^{ade}, Ulrika Axius^a, Werner J. Pichler^a, Daniel Yerly^a

^a ADR-AC GmbH, Adverse Drug Reactions – Analysis and Consulting, Bern, Switzerland

^b INRAAIC, Meyrin, Geneva, Switzerland

^c Clinical Immunology, Hôpital de la Tour, Meyrin, Switzerland

^d Löwenpraxis Luzern, Lucerne, Switzerland

^e Klinik St. Anna, Lucerne, Switzerland

	Antigen	Activation cut-off value (relative to negative control)	Convalescent (N=30)			Unexposed (N=10)			PPV	NPV
			+	-	Sensitivity	+	-	Specificity		
CD4 ⁺ T cells	S1	2	26	4	87%	1	9	90%	96%	69%
	S2							90%	96%	53%
	NC							100%	100%	63%
	Combined Statistics	2 in stimuli						90%	97%	82%
CD8 ⁺ T cells	Membrane Mix							70%	80%	28%
	Spike Mix							70%	83%	32%
	NC Mix							90%	94%	41%
	Combined Statistics	2 in stimuli condition	25	5	83%	4	6	60%	85%	46%

Combined
Sensitivity 93%
Specificity 90%

Analysis

Result

IMMUNOLOGY

SARS-CoV-2 Lymphocytes Analysis

T cells CD4

Spike 1

▲ positive

Spike 2

▲ positive

Nucleocapsid

negative

T cells CD8

Membrane mix

▲ positive

Spike mix

▲ positive

Nucleocapsid mix

▲ positive

B cells CD19

Spike 1

negative

Spike 2

negative

Nucleocapsid

negative

Observation

- Significant activation of CD4⁺ T cells detected in response to Spike 1 and Spike 2 antigens.
- Significant activation of CD8⁺ T cells in response to the tested peptide mixes of Membrane, Spike and Nucleocapsid.
- No activation of B cells.

Conclusion

Immunization to SARS-CoV-2 clearly detected.




Scan for details


Does vaccination induce T cell responses?

Helper T cells
CD4+



Vaccinated, N=9



Convalescent, N=4



Cytotoxic T cells
CD8+



✓ **Yes:** T cell responses are induced by vaccination

Clinical applications for Lymphocyte Analysis


1

- Persons with **immunodeficiencies/suppression**
 - Immunosuppression **treatments**
e.g. Transplantation, Cancer therapy, Biological treatments
 - Primary immunodeficiency

2

- Do T cells induced from vaccination & past infections still recognise **mutant** variants?
 - Omicron (2 years after original strain identified in Wuhan)

T cell immunity in **immunosuppressed** individuals



coloured points = antibody negative individuals

Healthy naïve controls N=10

Immunosuppressed (N=40)

- Vaccinated N=32
- Infected N=8

Even in immunosuppressed individuals:

- ✓ Vaccination induces both helper & cytotoxic T cell responses
- ✓ Past infection also has a trend towards positive helper and cytotoxic T cell responses
- ✓ but... Important to check at the **individual patient level** (*personalised* healthcare)



T cell immunity across **different immunosuppressants**

→ overall good vaccination-induced T cell responses

2° Immunodeficiency
N=20
after 2X vaccination

Secondary immunodeficiency			N=20
Target	Cells affected	Substances	N
CD20	B cells	Rituximab, Ocrevus, Obinutuzumab	7
Sphingosine 1-phosphate receptor (S1PR)		Fingolimod	4
Bruton's Tyrosine Kinase (BTK)		Ibrutinib	2
Calcineurin (→ IL-2)	T cells	FK506 (Tacrolimus)	1
CTLA-4		Belatacept	1
IL-12 & IL-23		Stelara (Ustekinumab)	1
General suppression	Lymphocytes	Imurek, Revlimid, Cellcept, Prednisolone	4

Healthy Donors
N=7
after 2X vaccination



Overall:

- ✓ Good T cell response
- ❖ but some poor responders

Recommendation:

- test for T cells
- extra caution / booster

Omicron – master of (antibody) evasion

Spike protein

(A) Delta






Image from: Kumar et al. (2022). *J. Med. Virol.* **94**(4)

(B) Omicron




mutation leading to
3D changes



antibody binding
impaired
→ viral evasion of
antibodies

degradation into linear peptides



T cells recognise linear peptides

✓ less affected by conformational changes
(recognition sites are conserved)

❖ when T cell epitopes are mutated,
problematic only for specific HLA-haplotypes

Stimuli adaptation for Omicron

Omicron mutations (60)




Image from: Chekol Abebe et al. (2022). *Infect. Drug Resist* 15

Wild-type Peptide sequence	Affected HLA alleles			Adapted for Omicron Peptide sequence
TEKSNIIRGW	B*44:02	B*44:03		IEKSNIIRGW
GEVFNATRF	B*40:01	B*44:02	B*44:03	DEVFNATRF
NSASFSTFK	A*23:01			NLAPFFETFK
KIADYNYKL	B*07:02			NIADYNYKL
AEHVNNNSY	B*44:03			AEYVNNNSY
SPRRARSVA	B*07:02			SHRRARSVA

- only 6 / 28 tested spike peptides affected
- most abundant HLA allele A*02 not affected

Helper
T cells

Cytotoxic
T cells

-	1	RPMI
-	2	RPMI + DMSO
+	3	Phytohaemagglutinin (PHA)
T	4	Spike 1 protein → Spike 1 RBD
T	5	Spike 2 protein
T	6	Nucleocapsid protein
T	7	Membrane peptides mix
T	8	Spike peptides mix
T	9	Nucleocapsid peptides mix

Is immunity still effective against Omicron?


Cohort

Vaccinated;

N=4


Infected (**from January 2022**) \pm Vax; N=6

likely Omicron infection
(+ ~80% Omicron)



T cells induced by vaccination / infection **still recognise Omicron mutant Spike protein**
 → likely **explains the reduced Covid-19 severity** despite increased breakthrough infections

Key Messages



Immunity exists as a
coordinated response

Immunity against SARS-CoV-2

Antibodies

Cellular immunity

Lymphocyte analysis

Detection of SARS-CoV-2 specific T lymphocytes possible

Test **available** and **in-use** for routine diagnostics

Application in routine diagnostics

- ✓ **Personalised** healthcare
- ✓ Immunodeficient patients:
Poor / No Antibody response –
→ **test for T cells**

Providing an unmet diagnostic tool

Confirmation of past-exposure / vaccination response

→ advice for caution when multiple immune parameters absent

→ preparedness for future mutants / other viruses

...any suggestions? Questions?

We are open to collaborate for further studies!

Acknowledgements



Dr.
Daniel Yerly



Dr. med.
Oliver Hausmann



Prof. Dr. med.
Werner Pichler



Dr. med.
Pierre Gumowski

ADR-AC team

Kevin Kammermann
Swelia Nussli
Benno Grabscheid
Florian Pichler

Sofia Blomquist
Zoë Jaccard

Blood donors